ENGLISH

Instructions for Use

DESCRIPTION

Avance® Nerve Graft is processed nerve allograft intended for bridging nerve discontinuities to support axon regeneration and to restore function. The cleansing process preserves the inherent microtubular structure of the native nerve extracellular matrix, whitening the cells, cellular debris and certain proteins such as chondroitin sulfate proteoglycans (CSPG).

Avance® Nerve Graft is implanted and connects the proximal and distal ends of a transected nerve. The processed nerve allograft serves as a scaffold allowing regenerating axons to grow into the patient's distal nerve tissue toward the target muscle or skin.

Avance® Nerve Graft offers structural characteristics and handling similar to nerve autograft, pliability of soft tissue, an intact epineurium to surround the nerve graft in place, and intact endoneurial tubes for the axons to grow through.

Avance® Nerve Graft is supplied sterile in a variety of lengths and diameters to allow the surgeon to match the size of the injured nerve and provide a tensionless repair of the defect. It is for single patient use only. Approximate graft lengths and diameters are listed on the package label.

REGULATORY CLASSIFICATION

Avance® Nerve Graft is processed and distributed in accordance with US FDA requirements for Human Cellular and Tissue-based Products (HCT/P) under 21 CFR Part 1271 regulations. US State regulations and the guidelines of the American Association of Tissue Banks (AATB). Additionally, international regulations are followed as appropriate.

Avance® Nerve Graft is to be dispensed only by or on the order of a licensed physician.

INDICATIONS FOR USE

Avance® Nerve Graft is processed nerve allograft (human) intended for the surgical repair of peripheral nerve discontinuities to support regeneration across the defect.

CONTRAINDICATIONS

Avance® Nerve Graft is contraindicated for use in any patient in whom soft tissue implants are contraindicated. This includes any pathology that would limit the blood supply and compromise healing or evidence of active infection.

WARNINGS

Careful donor screening, laboratory testing, tissue processing, and gamma irradiation have been utilized to minimize the risks of transmission of relevant communicable diseases to the patient. As with any processed human donor tissue, Avance® Nerve Graft cannot be guaranteed to be free of all pathogens.

Do not reuse or re-sterilize Avance® Nerve Graft and do not refreeze the graft once it has been thawed.

DONOR RECOVERY AND SCREENING

After consent for donation is obtained, surgical recovery of the peripheral nerve tissue is performed in an aseptic manner by appropriately licensed US State regulators. Donor eligibility is carefully evaluated as required by the US FDA and US State regulations. Additionally, donor eligibility is determined in accordance with AATB standards and appropriate international regulations. Tissue donors are evaluated for high risk behaviors and relevant communicable diseases. Evaluation includes a review of the donor medical and social history, a physical assessment of the donor at time of tissue recovery, a review of autopsy (if performed), serology testing, tissue recovery microbiology, and cause of death.

Each donor is tested and shown to be **negative or nonreactive** to the following:

- Human Immunodeficiency Virus (HIV) Type 1 Antibody
- Human Immunodeficiency Virus (HIV) Type 2 Antibody
- Hepatitis C Virus (HCV) Antibody
- Hepatitis B Virus (HBV) Surface Antigen
- Hepatitis B Virus (HBV) Core Antibody (total)
- Syphilis Rapid Plasma Reagin or Treponemal Specific Assay
- Human Immunodeficiency Virus (HIV) Nucleic Acid Test (NAT)
- Hepatitis C Virus (HCV) Nucleic Acid Test (NAT)

Additional testing may be performed, as required by local authorities in international markets, and may include the following:

- Hepatitis B Virus (HBV) Nucleic Acid Test (NAT)
- Human T-Cell Lymphotropic Virus (HTLV) Type 1 Antibody
- Human T-Cell Lymphotropic Virus (HTLV) Type 2 Antibody

All testing is performed by a laboratory certified to perform such testing on human specimens under the US Clinical Laboratory Improvement Amendments of 1988. The testing is conducted using test kits approved by the US FDA for testing cadaveric specimens where available.

The Medical Director of AxxoGen (US state licensed) has determined that the tissue is suitable for transplantation in humans. Records of all testing and medical releases are maintained by AxxoGen.

PROCESSING

Avance® Nerve Graft is processed in controlled environments using Good Tissue Practice (GTP) methods designed to prevent contamination and cross contamination of the tissue. Processing involves the use of proprietary physiological buffers, enzyme and surfactants and the processing tissue and processing processing agents. The cleansing process preserves the inherent microtubular structure of the native nerve extracellular matrix while clearing the cells, cellular debris and certain proteins such as chondroitin sulfate proteoglycans (CSPG). After completion of processing, Avance® Nerve Graft is sized, packaged and sterilized using gamma irradiation in accordance with ISO 11137 guidelines.

Avance® Nerve Graft has been tested in accordance with ISO 10993 standards. The test results demonstrated that the processed nerve allograft is biocompatible.

The processed nerve allograft tissue will naturally vary in color from white, off-white, pink, pale pink, and yellow to pale yellow. Occasional dark spots or localized discoloration is a normal occurrence.

HOW SUPPLIED

Avance® Nerve Graft is placed into a plastic tray and then inserted into chevon pouches. Each chevon pouch is heat-sealed to provide a sterile barrier and chevon pouch has a chevon seal. The outer chevon pouch is full to provide a moisture barrier. Approximate graft lengths and diameters are listed on the package label. Avance® Nerve Graft is irradiated and supplied frozen. Contents of the foil package are sterile unless the package is opened or damaged.

TRANSPORT AND STORAGE

Avance® Nerve Graft is shipped frozen on dry ice via validated transportation methods. Upon arrival, the product should be removed from the shipping container and placed in a freezer at or below -40°C (-40°F). Expiration date is three (3) years from date of packaging provided that the Avance® Nerve Graft has been stored at temperatures at or below -40°C (-40°F). See product label for expiration date.

Temporary storage conditions

Temporary storage of Avance® Nerve Graft at -20°C to -40°C (-4°F to -40°F) is limited to six (6) months total, and grafts stored at this temperature range must then be transferred to -40°C (-40°F) or colder freezer, used or discarded.

It is the responsibility of the health care institution to track the expiration date of the Avance® Nerve Graft and ensure that the product is stored properly.

WARNING

If the outer foil chevon pouch and/or inner Tyvek® chevon pouch is compromised (shows evidence of being torn or opened in any manner), **DO NOT USE** the Avance® Nerve Graft and notify AxxoGen Customer Care immediately.

INSTRUCTIONS FOR USE

- Follow standard operating procedures for exposure and mobilization of the injured peripheral nerve.
- Determine the injured nerve diameter in millimeters (mm) using a suitable measuring instrument.
- Select Avance® Nerve Graft(s) of comparable diameter to match the native nerve and of sufficient length to ensure a tensionless free repair.
- To prepare the Avance® Nerve Graft:
 - Remove the foil chevon pouch containing Avance® Nerve Graft. Instructions for Use, patient record labels and Tissue Utilization Report (TUR) from the package.
 - Compare the distinct lot number on the foil chevon pouch with the lot number on the package. If the numbers do not match, **DO NOT USE** the product and notify AxxoGen Customer Care immediately.
 - Using standard aseptic technique, peel open the outer foil chevon pouch and pass the inner Tyvek® chevon pouch to the sterile field for further handling.
 - Open the Tyvek® chevon pouch and remove the plastic tray.
 - Open the plastic tray and fill pre-molded thermal reservoir with room temperature sterile saline or sterile Lactated Ringers Solution (LRS). Do not heat the graft or add heated saline or LRS to the graft.
 - Allow Avance® Nerve Graft to thaw completely before use, which will take approximately 5 - 10 minutes. Once thawed, Avance® Nerve Graft is soft and pliable throughout. Avance® Nerve Graft must be either implanted or discarded within 12 hours. **NEVER IMPLANT A PARTIALLY OR FULLY FROZEN PRODUCT.**
- Handle Avance® Nerve Graft by outer most epineurium and avoid crimping or crushing the graft.
- Implant Avance® Nerve Graft using the same tensionless surgical technique when implanting a nerve autograft. Either end of the processed nerve allograft can be coapted to the proximal stump of the host nerve.
- Destroy any thawed allograft tissue not used in the surgical procedure in accordance with local, state and federal or country regulations for disposal of human tissue.
- Complete and send the Tissue Utilization Report (TUR) back to AxxoGen.

TISSUE UTILIZATION REPORT (TUR)

Each Avance® Nerve Graft package contains a Tissue Utilization Report (TUR). In accordance with US FDA, US Joint Commission and international requirements, a

TUR should be completed for each Avance® Nerve Graft used in the procedure and returned to AxxoGen or other representative as described on the TUR.

Record the distinct HCT/P identification code in hospital or facility records and in the patient's file. Complete all information on the cart, attach ONE (1) peel-off label of each Avance® Nerve Graft. Split and send return to AxxoGen or other representative as described on the TUR.

It is the responsibility of the health care institution to maintain recipient records for the purpose of tracking tissue post-implantation. The Tissue Utilization Report is a NOT intended to be a substitute for a facility's internal tissue transplantation tracking system.

POTENTIAL COMPLICATIONS

The following adverse events are anticipated and may potentially occur during treatment:

- mild incisional redness;
- tenderness of surgical areas;
- mild edema of surgical areas;
- controllable pain at surgical area;
- decreased sensation; and,
- numbness

These adverse events are considered expected and are not required to be recorded unless they increase in severity.

Inherent risks of any surgical procedure include, infection, blood loss, and anesthesia associated complications. Complications specific to nerve reconstruction include pain, decreased or increased sensitivity, and impaired motor or sensory function. As with all peripheral nerve surgery, there is a risk of failure of the nerve to regenerate.

Hypersensitivity, allergic reactions, or other adverse immune responses have not been observed in preclinical studies or reported clinically with the use of Avance® Nerve Graft. Because Avance® Nerve Graft is composed of proteins such as collagen and laminin, the potential may exist for such reactions.

Avance® Nerve Graft is processed human nerve tissue. As with all donated human tissue products the risk for transmission of communicable disease does exist. Robust donor screening and selection criteria, completed as required by the FDA and in accordance with AATB, state, and federal guidelines, provide control and assurance of sterilization with gamma irradiation greatly reduce but cannot totally eliminate this risk. As disease screening methods are limited, certain diseases may not be detected. The following complications of tissue transplantation may occur:

- Transmission of diseases of unknown etiology;
- Transmission from known infectious agents including, but not limited to viruses, bacteria, and fungi;

DISPOSAL

Dispose of Avance® Nerve Graft in accordance with local, state and federal or country regulations for disposal of human tissue.

REFERENCES

- ISO 11317:2006 Sterilization of health care products – Radiation guidelines
- ISO 10993:2003 Biological evaluation of medical devices

INQUIRIES

For additional information, to place an order, or to report errors, accidents or adverse reactions, contact:

AxxoGen Customer Care
Phone: 888-AxxoGen (888.296.4361) or 386.462.6800
Email: CustomerCare@AxxoGenInc.com

Customers outside of US: Contact the AxxoGen authorized distributor servicing your facility, email AxxoGen Customer Care or contact AxxoGen directly in the US at 386.462.6800.

Symbols used on Product Packaging

REF LOT

Product Code

Lot Number

Temperature Limit

Do Not Reuse/Restertize

Contents are sterile unless outer package is damaged, Radiation Sterilization

See Instructions For Use

Expiration Date

This product, its manufacture and use is covered by one or more of U.S. Patent No. 6,972,168, U.S. Patent No. 7,402,319, U.S. Patent No. 7,732,200, Japanese Patent No. 4749667, Europe Patent No. EP1425390 and other patents pending

Canada CTO Registration #100065

Avance® Nerve Graft is processed in the United States by:

AxxoGen

13631 Progress Blvd., Suite 400
Alachua, FL 32615
www.AxxoGenInc.com

Avance® Nerve Graft is a registered trademark of AxxoGen. Tyvek® is a registered trademark of E.I. DuPont De Nemours & Co.

Avance® Nerve Graft is a registered trademark of AxxoGen.

The processed nerve allograft tissue will naturally vary in color from white, off-white, pink, pale pink, and yellow to pale yellow. Occasional dark spots or localized discoloration is a normal occurrence.

HOW SUPPLIED

Avance® Nerve Graft is placed into a plastic tray and then inserted into chevon pouches. Each chevon pouch is heat-sealed to provide a sterile barrier and chevon pouch has a chevon seal. The outer chevon pouch is full to provide a moisture barrier. Approximate graft lengths and diameters are listed on the package label. Avance® Nerve Graft is irradiated and supplied frozen. Contents of the foil package are sterile unless the package is opened or damaged.

TRANSPORT AND STORAGE

Avance® Nerve Graft is shipped frozen on dry ice via validated transportation methods. Upon arrival, the product should be removed from the shipping container and placed in a freezer at or below -40°C (-40°F). Expiration date is three (3) years from date of packaging provided that the Avance® Nerve Graft has been stored at temperatures at or below -40°C (-40°F). See product label for expiration date.

Temporary storage conditions

Temporary storage of Avance® Nerve Graft at -20°C to -40°C (-4°F to -40°F) is limited to six (6) months total, and grafts stored at this temperature range must then be transferred to -40°C (-40°F) or colder freezer, used or discarded.

It is the responsibility of the health care institution to track the expiration date of the Avance® Nerve Graft and ensure that the product is stored properly.

WARNING

If the outer foil chevon pouch and/or inner Tyvek® chevon pouch is compromised (shows evidence of being torn or opened in any manner), **DO NOT USE** the Avance® Nerve Graft and notify AxxoGen Customer Care immediately.

INSTRUCTIONS FOR USE

- Follow standard operating procedures for exposure and mobilization of the injured peripheral nerve.
- Determine the injured nerve diameter in millimeters (mm) using a suitable measuring instrument.
- Select Avance® Nerve Graft(s) of comparable diameter to match the native nerve and of sufficient length to ensure a tensionless free repair.
- To prepare the Avance® Nerve Graft:
 - Remove the foil chevon pouch containing Avance® Nerve Graft. Instructions for Use, patient record labels and Tissue Utilization Report (TUR) from the package.
 - Compare the distinct lot number on the foil chevon pouch with the lot number on the package. If the numbers do not match, **DO NOT USE** the product and notify AxxoGen Customer Care immediately.
 - Using standard aseptic technique, peel open the outer foil chevon pouch and pass the inner Tyvek® chevon pouch to the sterile field for further handling.
 - Open the Tyvek® chevon pouch and remove the plastic tray.
 - Open the plastic tray and fill pre-molded thermal reservoir with room temperature sterile saline or sterile Lactated Ringers Solution (LRS). Do not heat the graft or add heated saline or LRS to the graft.
 - Allow Avance® Nerve Graft to thaw completely before use, which will take approximately 5 - 10 minutes. Once thawed, Avance® Nerve Graft is soft and pliable throughout. Avance® Nerve Graft must be either implanted or discarded within 12 hours. **NEVER IMPLANT A PARTIALLY OR FULLY FROZEN PRODUCT.**
- Handle Avance® Nerve Graft by outer most epineurium and avoid crimping or crushing the graft.
- Implant Avance® Nerve Graft using the same tensionless surgical technique when implanting a nerve autograft. Either end of the processed nerve allograft can be coapted to the proximal stump of the host nerve.
- Destroy any thawed allograft tissue not used in the surgical procedure in accordance with local, state and federal or country regulations for disposal of human tissue.
- Complete and send the Tissue Utilization Report (TUR) back to AxxoGen.

TISSUE UTILIZATION REPORT (TUR)

Each Avance® Nerve Graft package contains a Tissue Utilization Report (TUR). In accordance with US FDA, US Joint Commission and international requirements, a

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