

The background of the slide is a grayscale scanning electron micrograph (SEM) showing numerous small, angular, and polyhedral particles, characteristic of silicon nitride. The particles vary in size and are distributed across the field of view. In the top left corner, the text '7-20A' is visible. In the bottom right corner, a scale bar is labeled '1 μm'.

# An Overview of Silicon Nitride as a Novel Biomaterial

Bryan J. McEntire, MBA,  
Alan Lakshminarayanan, Ph.D.,  
B. Sonny Bal, MD, MBA, JD, and  
Thomas J. Webster, Ph.D.



A M E D I C A<sup>®</sup>

# $\text{Si}_3\text{N}_4$ has many industrial uses




## Bearings

- Machine tool spindles
- Dental handpieces
- Space shuttle/satellites
- Bicycles/in-line skates




## Cutting tools



## Turbo-machinery

- Turbine rotors
- Stators
- Glow plugs



## Electronics

# $\text{Si}_3\text{N}_4$ medical applications

## Spine



☐ Cages and corpectomy devices

## Reconstructive



☐ Total joint arthroplasty



A M E D I C A<sup>®</sup>



# Si<sub>3</sub>N<sub>4</sub> Advantages

---

- △ High flexural strength
- △ High fracture toughness
- △ Reasonable Weibull modulus
- △ Low wear rate
- △ Biocompatible
- △ Possesses antibacterial behavior
- △ Textured or highly polished surfaces
- △ Allows for bone on-growth
- △ Allows for bone in-growth
- △ Phase stable composition
  - Immune to hydrothermal degradation

# Si<sub>3</sub>N<sub>4</sub> Advantages

## Comparative Properties of Biomaterials<sup>1</sup>

Property	Units	Si <sub>3</sub> N <sub>4</sub>	Al <sub>2</sub> O <sub>3</sub>	ZTA	YSZ	CoCr	PEEK Optima®	Ti - Alloy	Cortical Bone
Density	g/cc	3.15-3.26	3.986	4.37	6.04	~8.5	1.29	4.43	1.85
Elastic Modulus	GPa	300-320	400-450	350	210	210-250	4.2	105-115	12-Aug
Compressive Strength	MPa	2500-3000	2000-5000	4300	2200	600-1800	130-140	950-990	130-190
Flexural Strength	MPa	800-1100	300-600	1000-1150	1050	-	160-180	-	-
Fracture Toughness	MPa·m <sup>1/2</sup>	8-11	4-5	5.7	10.5	50-100	-	75	-
Hardness	GPa	13-16	14-16	19.1	12.5	3-4	-	3.4	-
Surface Composition	NA	SiNH <sub>2</sub> & SiOH	Al <sub>2</sub> O <sub>3</sub>	ZrO <sub>2</sub> , Al <sub>2</sub> O <sub>3</sub>	ZrO <sub>2</sub>	CoO, Cr <sub>2</sub> O <sub>3</sub>	-OH Groups	TiO <sub>2</sub> , Al <sub>2</sub> O <sub>3</sub>	Ca <sub>5</sub> (PO <sub>4</sub> ) <sub>3</sub> (OH)
Surface Charge at pH =	NA	Positive	Positive	Positive	Positive	Negative	Negative	Negative	-



1. B. S. Bal, et al., "Orthopedic applications of silicon nitride ceramics," *Acta Biomaterialia*, 8, (2012), 2889-2898.

30KV


X10,000

6mm




# Bioceramic Microstructures


$\text{Al}_2\text{O}_3$  – Polished




ZTA – Polished



$\text{Si}_3\text{N}_4$  – Polished



$\text{Si}_3\text{N}_4$  – As-Fired



A M E D I C A<sup>®</sup>

# Crack Propagation in Ceramics

## Conventional ceramics

Low fracture toughness due to non-torturous crack path



## $\text{Si}_3\text{N}_4$ ceramics

High fracture toughness results from a torturous crack path



AMEDICA®

30KV X10,000

5mm

# Phase and Strength Stability of $\text{Si}_3\text{N}_4$



Autoclave conditions:  
120°C, 1 atm steam for  
up to 100 hrs

No change in  
phase composition

No change in  
flexural strength

Inherently stable microstructure


B. S. Bal, et al., "Fabrication and Testing of Silicon Nitride Bearings in Total Hip Arthroplasty," *J. Arthroplasty*, 24, [1], 110-116, (2009).





# Si<sub>3</sub>N<sub>4</sub> Femoral Head Burst Test Results

Strength of 28 mm **MC<sup>2</sup>** femoral heads in various designs



# Si<sub>3</sub>N<sub>4</sub> Biocompatibility


Si <sub>3</sub> N <sub>4</sub> Biocompatibility Tests (ISO 10993)	Pass
Cytotoxicity	✓
Systemic Toxicity	✓
Sensitization	✓
Genotoxicity	✓
Hemolysis	✓
Muscle Implantation (2 & 4 wk.)	✓
Sterilization compatibility (steam and γ - irradiation)	✓

## Confirmation

- Kue R, Sohrabi a, Nagle D, Frondoza C, Hungerford D., “Enhanced proliferation and osteocalcin production by human osteoblast-like MG63 cells on silicon nitride ceramic discs,” *Biomaterials*. 1999;20(13):1195–201.
- Neumann A, Reske T, Held M, et al., “Comparative investigation of the biocompatibility of various silicon nitride ceramic qualities in vitro,” *Journal of Materials Science. Materials in Medicine*. 2004;15(10):1135–40.

A M E D I C A<sup>®</sup>

# $\text{Si}_3\text{N}_4$ Wear Characteristics



28mm  
**PE/CoCr**  
#1-3

28mm  
**PE/SN**  
#1-3


28mm  
**XLPE/SN**  
#1-3

40mm  
**XLPE/SN**  
#1-3

B. J. McEntire, "Hip Simulator Wear Testing of Infinia® Femoral Heads,"  
Amedica Internal Research Report #RR 10001-11 A-2, (June 17, 2012).




# $\text{Si}_3\text{N}_4$ Wear Characteristics



B. S. Bal, et al., "Fabrication and Testing of Silicon Nitride Bearings in Total Hip Arthroplasty," *J. Arthroplasty*, 2009;24(1):110-116.

# $\text{Si}_3\text{N}_4$ Imaging Characteristics




J. Bernero, et al., "Medical Imaging Characteristics of Silicon Nitride," SAS Conference, Miami, (2008).



# Wettability of Various Biomaterials

## Contact Angle Measurements for Bioceramics<sup>1</sup>




## Contact Angle Measurements for Biomaterials<sup>2</sup>

Material	Water Contact Angle
PEEK	95°
Titanium	76°
$\text{Si}_3\text{N}_4$	39°

W = Water; BS = Diluted Bovine Serum;  
PHS2 = Hank's Balanced Salt Solution

1. M. Mazzocchi et al., "On the possibility of silicon nitride as a ceramic for structural orthopaedic implants. Part II: Chemical stability and wear resistance in body environment," *J. Mater. Sci. Mater. Med.*, (2008);19:2889-2901.
2. D. Gorth et al., "Decreased bacteria activity on  $\text{Si}_3\text{N}_4$  surfaces compared with PEEK or titanium," *Int. J. Nanomedicine*, (2012), in press.

# Porous and Dense $\text{Si}_3\text{N}_4$ Constructs




**C<sup>SC</sup>**®

**$\text{Si}_3\text{N}_4$  Cancellous Structured Ceramic**

**MC<sup>2</sup>**®

**$\text{Si}_3\text{N}_4$  Micro-Composite Ceramic**



- ▲ Strength
- ▲ Superior Imaging
- ▲ Osteo-integration / Conductive Scaffold
- ▲ Bacterial Resistant Material




A M E D I C A®

# Bone Grows into C<sup>s</sup>C<sup>®</sup> Si<sub>3</sub>N<sub>4</sub>

## Sheep Study – Medial femoral condyle

*(Cancellous Structured Ceramic (C<sup>s</sup>C<sup>®</sup>), ~72% Porous Si<sub>3</sub>N<sub>4</sub> Construct)*

Retrieved implant  
12 wks post-op



Bone penetrated >3mm  
into CSC<sup>®</sup>



M. C. Anderson, et al., "Bone In-Growth into Porous Silicon Nitride," *J. Biomed. Mat., Part A*, (2010);92A(4);1598-1605.



# Clinical Results of C<sup>s</sup>C<sup>®</sup> Si<sub>3</sub>N<sub>4</sub>

## △ Currently Used as an Implant in the EU

- One Level Cervical Fusion Without Instrumented Fixation
- Excellent Fusion and Osteointegration





Mark P. Arts, MD, PhD, Neurosurgeon, *Medical Center Haaglanden*, The Netherlands, Personal Communication, (May, 2012)

## MRI – Three Months Post-Operatively







# Texturing Leads to Osteointegration



- △ Micro-roughened surface
- △ Increases surface area for on-growth
- △ Inter-digitated bony fixation



A M E D I C A




# Anti-Infective Properties of Si<sub>3</sub>N<sub>4</sub>

## *In vitro* studies

### *Staph. epi* Biofilm Production<sup>1</sup>

(Five strains of bacteria examined on three biomaterials)



1. D. Gorth et al., "Decreased bacteria activity on Si<sub>3</sub>N<sub>4</sub> surfaces compared with PEEK or titanium," *Int. J. Nanomedicine*, (2012), in press.


A M E D I C A<sup>®</sup>

# Anti-Infective Properties of Si<sub>3</sub>N<sub>4</sub>

## *In vitro* studies

### *Staph. epi* Count of Live Bacteria<sup>1</sup>

(Si<sub>3</sub>N<sub>4</sub> demonstrates remarkable resistance to biofilm formation and growth)




1. D. Gorth et al., "Decreased bacteria activity on Si<sub>3</sub>N<sub>4</sub> surfaces compared with PEEK or titanium," *Int. J. Nanomedicine*, (2012), in press.


A M E D I C A<sup>®</sup>

# *In vivo* Wistar Rat Calvarial Study<sup>1</sup>


**PEEK** - 5% new bone and 95% bacteria at implant interface; 21% new bone growth and 88% bacteria growth in surgical area




**Si<sub>3</sub>N<sub>4</sub>** – 65% new bone at implant interface, 71% new bone in surgical area



**Ti** – 19% new bone at implant interface, 36% new bone in surgical area



**Si<sub>3</sub>N<sub>4</sub>** – 52% new bone at implant interface, 66% new bone in surgical area



## Histology

3 months post-operatively. no bacteria inoculation


1. T. J. Webster et al., "Anti-infective and osteointegration properties of silicon nitride, poly(ether ether ketone), and titanium implants, *Acta Biomaterialia*, (2012), in press.

A M E D I C A<sup>®</sup>




# *In vivo* Wistar Rat Calvarial Study<sup>1</sup>


**PEEK** - 8% new bone at implant interface, 24% new bone in surgical area




**Si<sub>3</sub>N<sub>4</sub>** - 25% new bone and 0% bacteria at implant interface, 39% new bone growth and 0% bacteria in surgical area



**Ti** - 9% bone and 67% bacteria at implant interface, 26% new bone and 21% bacteria in surgical area



**Si<sub>3</sub>N<sub>4</sub>** - 21% new bone and 0% bacteria at implant interface, 42% new bone and 0% bacteria in surgical area



## Histology


3 months post-operatively with 10<sup>4</sup> *S. epi.* bacteria inoculation

1. T. J. Webster et al., "Anti-infective and osteointegration properties of silicon nitride, poly(ether ether ketone), and titanium implants, *Acta Biomaterialia*, (2012), in press.

A M E D I C A<sup>®</sup>

# *In vivo* Wistar Rat Calvarial Study<sup>1</sup>

## Push-out strengths without bacteria inoculation




1. T. J. Webster et al., "Anti-infective and osteointegration properties of silicon nitride, poly(ether ether ketone), and titanium implants, *Acta Biomaterialia*, (2012), in press.

A M E D I C A<sup>®</sup>



# *In vivo* Wistar Rat Calvarial Study<sup>1</sup>

Push-out strengths **with**  $10^4$  *Staph. epi.* inoculation



1. T. J. Webster et al., "Anti-infective and osteointegration properties of silicon nitride, poly(ether ether ketone), and titanium implants, *Acta Biomaterialia*, (2012), in press.

A M E D I C A



# Summary

**Si<sub>3</sub>N<sub>4</sub>**

High flexural strength and phase stability

High fracture toughness; Resists cracking and stress risers

Chemically stable; Resistant to biodegradation

Biologically inert; biocompatible

Highly wear resistant; Low wear rates

Possesses optimal imaging qualities

Excellent wettability to biologic fluids

Good bone *in-growth* in the form of CSC<sup>®</sup>

Good bone *on-growth* in the form of MC<sup>2</sup><sup>®</sup>

Leads to bone apposition instead of a fibrous layer

Possesses inherent anti-infective characteristics

A M E D I C A<sup>®</sup>